## L790 Biological Microscope



**L790** series biological microscopes are equipped with monocular, achromatic objectives and wide field eyepieces. They have provides clear image and convenience operation. They can be used for medicine agriculture area widely. They are the ideal instruments in medical treatment, teaching demonstration.

## Features

- ▲ With wide field eyepieces and achromatic objectives, the field of view is clear.
- **Δ** Coaxial coarse/fine focus system. Minimum division of fine focusing: 4μm.
- ▲ High brightness white LED lamp transmission illumination

### **Standard Configuration**

Model	L790A	L790B
Eyepiece	WF10X(Φ18mm)	,
	Achromatic 4X/0.10	Achromatic 4X/0.10
Objective	Achromatic 10X/0.25	Achromatic 10X/0.25
	Achromatic 40X/0.65(Spring)	Achromatic 40X/0.65(Spring)


Eyepiece tube	Achromatic100X/1.25(Spring oil)  Monocular,360° rotatable (Inclination of 45°)		
Nosepiece	Triple(Frontward ball bearing outer locating)	Quadruple (Frontward ball bearing outer locating)	
Focus system	Coaxial coarse/fine focus system. Minimum division of fine focusing: 4µm		
Stage	Fixed(Size:120mmX120mm, With five holes fixed iris diaphragm)	Double-layer mechanical, size:125mmX115mm, move range:70mmX30mm	
Abbe condenser		Rack &pinion adjustment	
Collector	For LED lamp		
Light source	LED High brightness white LED lamp		

## **Optional accessories**


Name	Sort/Technique parameter	Sort/Technique parameter	
	Wide field WF16X(Φ11m	Wide field WF16X(Φ11mm)	
Eyepiece	Point 10X(Φ18mm)	X(Φ18mm)	
	Fixed(Size:120mmX120mm	Fixed(Size:120mmX120mm, With single lens and iris diaphragm)	
Stage	Fixed(Size:120mmX120mm, with slide moving holder and Abbe condenser)		040004
Filter	Blue filter	Blue filter	
	Green filter	Green filter	
	Yellow filter	Yellow filter	
	Frosted filter	Frosted filter	
	For Fluorescent and Haloge	For Fluorescent and Halogen lamp	
Collector	For Incandescent lamp	gen lamp	060004
Light source	High brightness white LED	With adapter	073011
	High brightness white LED	With adapter and charge batteries	073012
	Halogen lamp	12V 20W	070019
	Fluorescent lamp	110V/220V, 5W Can be choose	071006
	Incandescent lamp	220V/110V, 20W Can be choose	072004


Camera	Photography	USB Output VIDEO Output	1096001 1108001
Adapter	Input: 230V;«50/60HZ (	Output: 9V 500mA	101001
(For LED illumination)	Input: 110V; «50/60HZ (	Output: 9V 500mA	101002

## Diagram



#### **Dimensions**





Items and specification subject to change, please contact ALLTION for current product information

### MICROSCOPES ARE OUR BUSINESS & SERVICE IS OUR PRIDE

# ALLTION (Wuzhou) Co., Ltd.

Address: Alltion Building, NO. 10, 3<sup>rd</sup> Road Wuzhou Industrial Park, Wuzhou, Guangxi, China 543000

Tel: +86-774-2836101 Fax: +86-774-2836192 E-mail: sales@alltion.com or sales@alltion-microscope.com

Web site: <a href="http://www.alltion.com">http://www.alltion-microscope.com</a>